








SERRES


EN	SERRES SUCTION BAG SYSTEM MANUAL.....	4
DE	ANLEITUNG FÜR SERRES ABSAUGBEUTEL-SYSTEM.....	6
FR	MANUEL DU SYSTÈME DE POCHE D'ASPIRATION SERRES.....	8
ES	MANUAL DEL SISTEMA DE BOLSA DE ASPIRACIÓN SERRES.....	10
IT	ISTRUZIONI D'USO PER IL SISTEMA CON SACCHE DI RACCOLTA PER LIQUIDI ASPIRATI SERRES.....	12
SV	BRUKSANVISNING SERRES SUGPÅSESYSYSTEM.....	14
DA	BRUGERVEJLEDNING FOR SERRES SUGEPOSESYSTEM.....	16
NO	VEILEDNING FOR SERRES SUGEPOSESYSTEM.....	18
FI	SERRES IMUPUSSIJÄRJESTELMÄN KÄYTTÖOHJE.....	20
NL	HANDLEIDING SERRES-OPVANGZAKSYSTEEM.....	22
CS	PŘÍRUČKA PRO SYSTÉM ODSÁVACÍHO VAKU SERRES.....	24
SK	NÁVOD K SYSTÉMU ODSÁVACIEHO VAKU SERRES.....	26
TR	SERRES EMİŞ TORBASI SİSTEMİ KILAVUZU.....	28
HU	SERRES SZÍVÓZSÁKRENDSZER TÁJÉKOZTATÓ.....	30
PL	INSTRUKCJA OBSŁUGI SYSTEMU WORKÓW DO ODSYSANIA SERRES.....	32
PT	MANUAL DO SISTEMA COM SACOS DE ASPIRAÇÃO SERRES.....	34
SL	NAVODILA ZA UPORABO ASPIRACIJSKEGA SISTEMA SERRES.....	36
EL	ΕΓΧΕΙΡΙΔΙΟ ΣΥΣΤΗΜΑΤΟΣ ΣΑΚΟΥ ΑΝΑΠΠΟΦΗΣΗΣ SERRES.....	38
TH	คู่มือระบบถุงดูดของเหลว SERRES.....	40
LT	„SERRES“ SIURBIMO MAIŠELIŲ SISTEMOS VADOVAS.....	42


3.4


3.5


4.1


4.2


5


6.1


6.2

6.3

6.4

EN SERRES SUCTION BAG SYSTEM MANUAL

The Serres Suction Bag System is used during medical procedures to collect liquids and secretions for disposal from patients. **Please read this manual carefully before installing the product and beginning the procedure.** The Serres suction bag system may only be installed and used by persons familiar with the product and used for the purpose stated by the manufacturer. Usage contrary to instructions is forbidden. The manufacturer cannot be held liable if the product is used in a manner and/or for purposes other than those described in these instructions. The system is not approved for use in the procedures where suctioned liquid is injected back to the patient. The Suction Bag is not approved for sample or specimen collection purposes but only for collection and disposal. Before using the product, always ensure that all of the system's parts are intact. Use of faulty products is strictly forbidden.

Parts of the system

The basic parts of the system include the suction bag **Figure A** (disposable), the suction canister and the grey, right-angled connector of the canister **Figure B** (re-usable). The volume of the suction bags and canisters is either 1,000 ml, 2,000 ml or 3,000 ml, depending on the product model. Parts of the suction bag system have been introduced in more detail in **Figure C**.

- C1** Angled connector (grey, attached to the suction canister)
- C2** Patient connector (white, attached to the suction bag)
- C3** Serial port (used when bags are connected in series)
- C4** Patient connection plug
- C5** Overflow protection
- C6** Lifting handle

There are different versions of the suction bag and their installation can differ from each other. Ensure that you have read the instructions related to the correct version.

WARNING! ONLY ORIGINAL SERRES PARTS AND ACCESSORIES MAY BE USED IN THE SUCTION BAG SYSTEM. **WARNING!** ALL SUCTION BAGS AND PRODUCTS THAT HAVE BEEN MARKED WITH THE  SYMBOL ARE DISPOSABLE AND RE-USING THEM IS STRICTLY FORBIDDEN.

Installation of the system

1. Installing the suction canister

Place the suction canister in a bracket in an upright position **Figure 1**. There are a range of brackets available for installing the suction canister (rail, wall, table and bed models). The suction canister can also be mounted to the canister holder in the movable Serres trolley. Connect the tubing from the vacuum source to the grey, angled connector at the rear of the canister.

2. Installing a single suction bag

2.1 Versions which have not been wrapped and taped. Unfold the suction bag and place it into the suction canister. **Figure 2.1**

2.2 Versions which have been wrapped and taped. Place the bag as it is in the suction canister **Figure 2.2** or follow the instructions in section 2.1.

WARNING! THE SUCTION BAG HAS TO BE INSTALLED INTO A SUCTION CANISTER OF EQUAL SIZE.

WARNING! ENSURE THAT THE SUCTION BAG FILM DOES NOT BECOME LODGED BETWEEN THE CANISTER AND THE LID.

3. Installing several suction bags (connection in series)

When large volumes of liquid are suctioned, Serres suction bags can be connected in series with serial tubes, vacuum tubes and T-connectors. **Figure 3.1**.

- Place the bags to suction canisters as in section 2 and attach the bags to the canisters by pressing the middle of the lid. **Figure 3.2**. If the Dual Filter Suction bag is used it must be placed to the last canister (3).

- Remove the grey, angled connector (canisters 2 and 3) and open the suction bags' serial ports (bags 1 and 2) as shown in **Figure 3.3**.


- Connect the suction canisters to each other with separate T-connectors and tubes in accordance with **Figure 3.4**. Use scissors to cut the tube to a suitable length.

NOTE! THE T-CONNECTOR AND VACUUM TUBE ARE RE-USABLE AND DO NOT NEED TO BE REPLACED BETWEEN OPERATIONS.

- Connect the serial tubes carefully through the opened serial port opening in the suction bag to the next suction bag's patient connector. **Figure 3.5**. For best performance, you are recommended to use a vacuum shift for the serial connections of the largest amounts of liquid.

4. Unfolding the suction bag

Suction bag is installed using vacuum. Turn the vacuum source on (**Figure 4.1**) and, at the same time, press the middle of the lid lightly. After the suction bag has straightened, close the patient connector using your finger so that the lid tightens to the suction canister. (**Figure 4.2**)

WARNING! BEFORE USE, ENSURE THAT THE VACUUM HAS BEEN CREATED AND THE SUCTION BAG IS FULLY INFLATED. REPEAT THE INSTALLATION IF THE FILM HAS NOT STRAIGHTENED.

Connect the patient tube to the patient connector. The system is ready for use.

The vacuum source can be turned off if the procedure is not immediately started.

WARNING! THERE IS A RISK OF CRACKING OF THE 1 L LID IF THE SUCTION BAG IS LEFT UNDER HIGH VACUUM FOR A LONG TIME. IF USING THE 1 L SUCTION BAG THE VACUUM SOURCE IS RECOMMENDED TO BE TURNED OFF IF THE PROCEDURE IS NOT SHORTLY STARTED. CHECK THAT THE PRODUCT IS INTACT BEFORE STARTING SUCTION. FAILURE TO DO SO CAN DELAY STARTING OF THE CRITICAL OPERATION.

Using the suction bag system

When the suction is on, the liquid flows inside the suction bag. The suction bag can be used throughout the procedure or until the suction bag is full. When the suction bag is full, overflow protection cuts the suction and prevents the flow of liquid to the suction source. The device is intended to measure quantity of liquids and secretions removed from the human body. The scale on the suction canister displays the volume of the suctioned substance in milliliters. The accuracy of the scale is +/- 100ml (only when the vacuum is on).

Persons using the Serres Suction Bag System should take appropriate safety precautions to protect themselves from coming into contact with patient liquids or secretions.

Do not attempt to move the trolley by pulling the patient tube.

Use of suction bags equipped with solidifying agent (Pre-gelled)

Some suction bag versions include a solidifying agent **Figure 5** that solidifies the suctioned liquid. Use pre-gelled suction bags similar to regular suction bags. When measuring the volume of the suctioned liquid, take note of the fact that the scale on the suction canister includes the amount of the solidifying agent contained in the suction bag.

1l pre-gelled suction bag: 25 ml

2l pre-gelled suction bag: 50 ml

3l pre-gelled suction bag: 85 ml

In addition, note the normal tolerance of the suction canister scale. If you use separate sachets or a solidifying agent sold as separate powder, use 25 g / 35 ml of substance per one litre of liquid. Add solidifier into the suction bag through the serial port and then close the port. Note the warnings in the solidifying agent packet. Gently squeezing the bag will help to mix the solidifier and the liquids promoting solidification.

NOTE! IF THE PRE-GELLED SUCTION BAGS ARE USED IN SERIAL CONNECTIONS, THE VACUUM SOURCE MUST BE ON FOR THE DURATION OF THE WHOLE PROCEDURE.
NOTE! SUCTION BAGS CONTAINING SOLIDIFYING AGENT MUST NOT BE EMPTIED INTO A SEWER OR IN SERRES NEMO EMPTYING DEVICE.

WARNING! DO NOT EAT SOLIDIFIER. AVOID INHALATION AND EYE CONTACT. SOLIDIFYING AGENT IS SLIPPERY WHEN WET AND MIGHT CAUSE SLIP HAZARD IF ACCIDENTALLY SPILLED ON THE FLOOR.

After the procedure

After the suctioning procedure, disconnect the patient tube and the patient connector (or alternatively straight connector or twin connector) and serial tube, if used **Figure 6.1**. Close the connection with the patient connector plug provided on the lid of the suction bag. In serial connections, also plug the serial port by pressing the center of the cap. **Figure 6.2**. Finally, turn off the vacuum source. **Figure 6.3**. Lift the suction bag from the handle **Figure 6.4**.

WARNING! DO NOT TURN OFF THE VACUUM SOURCE BEFORE YOU HAVE CLOSED THE SUCTION BAG.

NOTE! DO NOT DISCARD OR UNNECESSARILY REMOVE THE REUSABLE SUCTION CANISTER, GRAY ANGLED CONNECTOR OR SILICONE VACUUM TUBING.

Disposal of used products

Serres suction bags, collection cups, measuring cups and serial tubes are disposable and must be replaced patient-specifically. If the same patient undergoes long-term treatment, it is recommended that the suction bag be replaced at least every 24 hours. Reusable products must be decommissioned if they have become damaged or

no longer meet the manufacturer's intended operational specifications. Used or decommissioned products must be disposed of in accordance with recycling labelling and hospital-specific instructions. Suction bags containing solidifying agent must not be emptied into a sewer.

WARNING! REUSE OF DISPOSABLE PRODUCTS IS STRICTLY FORBIDDEN. REUSE REDUCES THE PERFORMANCE OF THE PRODUCT AND CAN CAUSE AN INFECTION RISK. WHEN HANDLING A USED SUCTION BAG, ATTENTION SHOULD BE PAID TO THE FACT THAT IT MAY CONTAIN INFECTIOUS WASTE.

Cleaning of re-usable products

The suction canister and angle connector can be washed (95°C) and autoclaved (121 °C). Remove the grey, angled connector before washing or autoclaving. Before autoclaving, make sure that the product does not contain any detergent remnants. Do not autoclave the product more than 30 cycles. Adapters and brackets (excluding 57815) can be washed (95°C). Other reusable products can be wiped with disinfectant.

Storage

Protect the packages from humidity, dirt and dust. Disposable products can be used for 5 years after the date on the label, excluding pre-gelled suction bags and solidifying agents that can be used for 2 years after the date on the label.

Reporting of serious incidents

Any serious incident related to the use of this product, should be reported to both the manufacturer and the health authority/competent authority where the product is installed.

Symbols used

LOT Batch code

REF Catalogue number

 Do not re-use

 Keep away from rain

 Keep away from sunlight

 Operating instruction

MD Medical Device

MFGDT (Manufacturing date)

 CE-mark

 DEHP-Free symbol

 Latexfree

 Quantity

 Refer to instruction manual booklet

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Serres Oy